



PNet control software for LAMBDA laboratory pumps, powder doser and gas flow controller


The PNet control software allows to realize flow profiles with the [LAMBDA laboratory peristaltic pumps PRECIFLOW, MULTIFLOW, HIFLOW, MAXIFLOW](#), the [syringe pump LAMBDA VIT-FIT](#), and also with the [LAMBDA DOSER powder dosing instrument](#) and the [LAMBDA MASSFLOW gas flow controller](#).

The instruments are connected to the computer through a RS-232 or RS-485 interface. Then, the COM port to which the instrument is connected and the address of the instrument are selected. Up to 6 LAMBDA instruments can be connected and controlled simultaneously. The PNet software runs under Windows NT, 2000, XP, Vista, 7, 8 and 8.1.

Main window

In this window you can:

- Enter set points for each instrument
- Configure, start and stop the set point profile
- Save the values of the set points in a file






Curves window (graphs)

This window displays the set point values of the different instruments as a function of time.



Profile Configuration (“Opt”-Button)

Here, set point profiles (up to 50 steps) can be introduced, saved and printed.



The 'Define a profile' dialog box contains the following table:

	Begin	End	Time (HHH:MM:SS)
1	999	0	000:09:59
2	999	150	001:00:00
3	150	200	001:00:00
4	0	0	000:00:00
5	0	0	000:00:00
6	0	0	000:00:00
7	0	0	000:00:00
8	0	0	000:00:00
9	0	0	000:00:00
10	0	0	000:00:00
11	0	0	000:00:00

Buttons: OK, Cancel, Open..., Save..., Print